

SafeSense

Digital Discovery Workbook

November 9th, 2018

Table of Contents

Value Proposition	3
FDA Regulatory Assessment	4
Use Case Definition	6
User Comps	8

Value Proposition

[Insert Here]

FDA Regulatory Assessment Questions

1. Who is the user? (e.g., Physician, nurse, patient)

Athletic Coach

Athletic Trainer

2. What is the purpose of the app? (Provide a 1 paragraph description)

The app will be used to display intensity mapping of impacts sustained by football players. The app will keep a record of lifetime impacts and provide real-time alerts if the impacts cross a pre-determined threshold. The pre-determined thresholds will be input based on information from published literature.

The alerts will provide a message stating a particular player(s) has crossed the pre-determined threshold. The athletic coach or trainer will then have the option to navigate the app and check the impacts sustained by the particular player(s). The display options would be:

- A graph showing the impacts (with intensity) sustained by the player(s) over a period of time
- A 3D image of the head showing impact intensities sustained by the player over a period of time

3. What are the inputs to the app? Be specific.

Date, atmospheric temperature, player statistics including name, height, weight, game position

4. What are the outputs of the app? Be specific.

Real time impact intensity mapping

5. If the app provides recommendation for patient treatment, briefly describe the scientific basis for the algorithm used to generate the treatment recommendations.

The app does not provide any recommendation for patient treatment.

6. What is the hardware platform?

App: ios/android

7. How will the software be provided? (e.g, downloaded from App store, pre-installed on hardware platform).

Downloaded from the app store.

Initial Regulatory Review

Donna-Bea Tillman, Ph.D, Biologics Consulting
December 8, 2018

FDA does not regulate protective gear such as football helmets or pads, even though the purpose of such equipment is to prevent injury. Therefore, I believe that if the app were only to provide “intensity mapping of impacts sustained by football players and maintain a record of lifetime impacts”, it would not meet the definition of a medical device, and thus would not be regulated by FDA.

The app is also intended to provide a message stating a particular player(s) has “crossed the pre-determined threshold” that will apparently be based on data from the literature. I presume that this threshold reflects the degree of impact that is likely to result in an injury to the player. This feature is a little less clear from a regulatory perspective. However, I think that it is probable that as long as the alert **does not specifically identify a particular disease or condition** (e.g., concussion or TBI), and instead simply informs the user that a potentially harmful impact has occurred, that the app would still probably not meet the definition of a medical device.

Note that FDA has classified software intended to be used to diagnose the presence of a concussion as a medical device, so it would be important to avoid making any direct reference to concussion or TBI.

In regard to the display options, as long as the 3-D image of the head is simply a general model and not an actual radiological image for the player in question, I think that displaying the impact intensities on the image would also not result in the app being regulated as a medical device.

Use Case Definition

Use case definition contains requirements for the SafeSense app.

Use Case Definition

Project: Safe Sense – Use Cases
Date: July 18, 2018
Owners: April Kwon, Sarah Jomaa

Goal:

1. *As an athletic trainer, I want to have a tablet application that comprehensively monitors and provides information, in real-time, on impacts football players experience on the field through helmet devices.*
2. *As a medical staff member, I want to have medical-related information collected from the helmet devices so I can make an accurate diagnosis more rapidly in an emergency.*

Coach Use Cases:

1. As a *trainer*, I want to see the home dashboard to show a summary of information about the individual players, the aggregated team data, and the devices that are connected.
2. As a *trainer*, before a game, I want to create a game view easily on the tablet with game participants' information and check the overall connection with devices.
3. As a *trainer*, during the game, I want to get a notice when a player gets a certain level of hits so I can decide whether I will take emergency measures or not.
4. As a *trainer*, I want a few different real-time graphs that compare different variables, such as the location of the hit or intensity level on each sensor on individual helmet gear, that way I can choose a view that I can check important changes on helmets.
5. As a *trainer*, I want to sort or filter collected real-time data for a different player, a group of players based on positions and to have different ranges of presentations of data so I can see the patterns based on different demographic parameters.
6. As a *trainer*, I want to see summary device-related data at the end of the game so I can see if there any need of additional treatment for a player.
7. As a *trainer*, I want to archive the sensor records during previous games so I can look at a player's data over the course of a play versus a season.
8. As a *trainer*, I want to change settings on which information or visualization to be shown in the dashboard.
9. As a *trainer*, I want to check the connection between and real-time data on the bluetooth connection and the battery status of the device so I can in a case replacement is needed.
10. As a *trainer*, I want to have a coach-medical team direct line on the application so I can be connected with the team more quickly in the emergency situation.

Medical staff Use Cases:

1. As medical staff, in an emergency situation, I want to be provided the location, the intensity, and the frequency of the hit(s) as quickly as possible so I can diagnose any possible concussion and take instant medical action.
2. As medical staff, I want to be informed about personal medical information, treatment history of a player that I should treat.


User Comps

User Comps are provided to allow an inventor the opportunity to review the navigation of the application's flow.

Low-Fi Prototype 1

Dashboard

Dashboard



Team Page




- 2 Click on each player to view their player profile.
- 1 Press and hold to add or remove players
- 4 Hold and drag over one half of the screen to view in split screen




View Ethmoid, Sphenoid, Occipital, Frontal, and Temporal sections of the skull. Tap once to view the slider and play/pause button, use them to go back and view a certain time.


Player Page



Games Page



Settings




Low-Fi Prototype 2

1) Home Dashboard (Coach/Trainer & Parent view)

1) Home Dashboard - Coach/Trainer View




2) Home Dashboard - Parent View




The home screen can be customized by users:
- Coach/trainer view: The summary of the team/selected players
- Parent view: A single player and related game information

2) Player Page (Add a player)




2) Player tab - Add a player



 Team information is stored in local database

 Contact information can connect to a player's coaches/physicians/parents


 The device can be connected to the application either during player profiling, or before a game

3) Game Page (Start game tracking)


3) Games - Add a Game



3) Games - Add a Game



3-1) During a game



The sign of the hit can be simplified




Real-time device activities (hit detection, device disconnection, etc.) are notified on the top menu




The filters can be based on the available parameters in player profiling page as well as the number/intensity of the hit detected.

1-2) Hit Visualization (Individual View)



4) Analytics & Setting



2) Setting

